



SIGMA Compact CsI Scintillation Spectrometer

Features

- All-in-one detector, electronics, and MCA
- Includes CsI scintillator detector with SiPM (No photomultiplier tube)
- Resolution better than 7.2% (at 662 keV)
- Small overall size, but big performance
 - SIGMA25: 35 x 35 x 105 mm; 200 g
 - SIGMA50: 35 x 35 x 130 mm; 300 g
- USB connection to PC
- USB powered, 250 mW max
- Interface to Genie™ 2000 software spectral analysis


Description

The SIGMA family is based on a CsI scintillator crystal coupled to a “state-of-the-art” SiPM (Silicon Photomultiplier) detector, high voltage, electronics and integrated MCA. All communication and power are provided over the mini-USB connector.

The SIGMA family comes in two scintillator sizes. The SIGMA25 has a 1 in³ CsI detector and the SIGMA50 has a 2 in³ CsI detector. All other aspects of the two devices are identical. Both units have a FWHM resolution specification of 7.2% at 662 keV, which is comparable to many NaI based solutions.



SIGMA integrated scintillator detector



Model	Case Size (mm)	Crystal Size (mm)	Volume cm ³	Resolution (%@Cs137)	Weight (g)
SIGMA50	35 x 35 x 130	25.4 x 25.4 x 51	32.8	<7.2	300
SIGMA25	35 x 35 x 105	25.4 x 25.4 x 25.4	16.4	<7.2	200

USB powered 250 mW

The profile of the SIGMA detectors makes them ideal for many portable spectroscopy applications where size, weight and power constraints are very important. The embedded MCA has a fixed number of channels, fixed high voltage and gain that provides a simple interface that is easy to setup and operate with essentially no configuration. There are three software options available for use with the SIGMA detector series, two of which are offered as free downloads from the CANBERRA website. All three packages are described on the next page.

SIGMA Compact CsI Scintillation Spectrometer

ISOCS™/LabSOCS™ Support


Quantitative measurements can be performed without the use of calibration sources using CANBERRA's mathematical efficiency calibration software ISOCS/LabSOCS. The SIGMA and GR1 family supports ISOCS efficiency calibrations. This characterization permits the user to calculate mathematical efficiencies using either ISOCS or LabSOCS software (sold separately). A unique characterization can be purchased for each specific detector for optimum accuracy or a generic characterization is available with slightly degraded accuracy specification.

Software Options

There are three software options available to communicate with the SIGMA family of sensors.

K-Spect Software


The K-Spect software can be downloaded free of charge from either the CANBERRA or Kromek websites. It establishes communications with the SIGMA enabling MCA control functions, spectrum acquisition, display, basic analysis and storage functions. Together, the SIGMA and K-Spect software establish a basic gamma-ray spectroscopy system that is suitable for qualitative spectral inspection and ROI based gamma analysis.



Genie 2000 Analysis

(through Kromek Connection Tool)

Although Genie 2000 cannot directly control the SIGMA family of sensors, CANBERRA is providing a simple control application that can be used to collect data and save it to a .CNF format file. The file can be analyzed using the full suite of Genie 2000 tools such as S501 Advanced Analysis and the ISOCS/LabSOCS mathematical efficiency packages. The Kromek Connection Tool is available as a free download from the CANBERRA website and is compatible with all current Genie 2000 analysis versions (Genie 2000 suite sold separately). The SIGMA, coupled with Genie 2000 analysis is the perfect combination for many portable spectroscopy applications requiring quantitative analysis (activity determination).



Kromek SDK

The SIGMA family of sensors is also supported by a Kromek software development kit (SDK). These software tools can be purchased from CANBERRA to allow end users to develop customized software applications to interact directly with the SIGMA. The tools provide complete setup and control of the key MCA functional parameters. The SDK supports both Windows® and Linux-based systems. Some programming knowledge is required to implement this product.

SIGMA Compact CsI Scintillation Spectrometer

Specifications

PERFORMANCE

- DETECTOR
 - SIGMA25 1" x 1" x 1" CsI(Tl) detector.
 - SIGMA50 1" x 1" x 2" CsI(Tl) detector.
- ENERGY RANGE – 50 keV to 1.5 MeV.
- ENERGY RESOLUTION – <7.2% FWHM at 662 keV (Standard Models).
- NUMBER OF CHANNELS – (fixed) 4096.
- MAXIMUM THROUGHPUT (USB) – 5000 cps.
- DIFFERENTIAL NON-LINEARITY – <±1%.

PHYSICAL

- SIZE
 - SIGMA25 – 34.5 x 34.5 x 105 mm.
 - SIGMA50 – 34.5 x 34.5 x 130 mm.
- WEIGHT
 - SIGMA25 – 200 g.
 - SIGMA50 – 300 g.

INPUTS/OUTPUTS

- USB (all models) – Signal and Power (Power consumption: 250 mW).

ENVIRONMENTAL

- OPERATING TEMPERATURE – 0 to 40 °C (32 to 104 °F).

ORDERING INFORMATION

SIGMA25	1" x 1" x 1" CsI (Tl) Spectrometer with – 7.2% FWHM resolution at 662 keV
SIGMA50	1" x 1" x 2" CsI (Tl) Spectrometer with – 7.2% FWHM resolution at 662 keV

Other Related Products

GR1

- GR1 Gamma-Ray Spectrometer with
- 2.5% FWHM resolution at 662 keV
 - Without MCX I/O ports

GR1+

- GR1+ Gamma-Ray Spectrometer with
- 2.0% FWHM resolution at 662 keV
 - Without MCX I/O ports

GR1-A

- GR1-A Advanced Gamma-Ray Spectrometer with
- 2.5% FWHM resolution at 662 keV
 - With MCX I/O ports (includes 3x MCX to BNC adaptors)

GR1-A+

- GR1-A+ Advanced Gamma-Ray Spectrometer with
- 2.0% FWHM resolution at 662 keV
 - With MCX I/O ports (includes 3x MCX to BNC adaptors)

ISXCZT-GR1

- Generic ISOCS Characterization for the GR1, GR1+, GR1-A or GR1-A+

ISXCSI25

- Generic ISOCS Characterization for the SIGMA25

ISXCSI50

- Generic ISOCS Characterization for the SIGMA50

LAB CPG (10 x 10 x 10)

- 1cm³ CZT detector for use with external MCA electronics

ISXCZT-LAB1000

- Generic ISOCS Characterization for the LAB CPG (10 x 10 x 10)

TN-15

- USB enabled Thermal Neutron Detector



Genie, ISOCS and LabSOCS are trademarks and/or registered trademarks of Mirion Technologies, Inc. and/or its affiliates in the United States and/or other countries.

All other trademarks are the property of their respective owners.

©2017 Mirion Technologies (Canberra), Inc. All rights reserved.

Copyright ©2017 Mirion Technologies, Inc. or its affiliates. All rights reserved. Mirion, the Mirion logo, and other trade names of Mirion products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners.

CANBERRA