



Nal Detector Accessories

Model 7419 Shield/Collimator



- All attitude
- Compact – swivel legs and folding handle
- Removable collimator – 44 mm (1.75 in.)
- Removable tin filter – 1 mm thick (0.040 in.)
- Lead shielding – 8 mm wall x 122 mm long (0.32 x 4.8 in.)
7419BS – 8 mm wall x 76 mm long (0.32 x 3.0 in.)
- Shipping weight – 6.5 kg (14 lb)
- Five models – specify when ordering
 - 7419A for 802-2 x .5
 - 7419B for NAID
 - 7419BS for NAID (with detachable legs)
 - 7419C for 802-2 x 2
 - 7419E for IPROS-2

