





Model INSITPOD Lightweight Collapsible Tripod




- Detector height adjustable ~35 – ~100 cm.
- Weighs 3.2 kg (7 lb).
- Legs stay in fixed position for quick and easy movement between measurement points.
- Small shipping package:
 - Legs easily removed and collapse to fit in 60 cm long package.
 - Attachment plate can stay connected to detector carrier.
- Accommodates both MAC and Big MAC.
- Accommodates both Slimline or Flanged preamplifiers, with or without RDC.
- Will not carry any additional weight (e.g. shield/collimators).

Model 717 Portable Detector Shield



- For MAC or Big MAC.
- 89 mm I.D. x 165 mm long (3.5 x 6.5 in.).
- 25 mm (1.0 in.) thick lead.
- Brass lined-steel jacketed.
- Shipping weight - 34.5 kg (76 lb).

Model 7411 Detector Submarine



- Accommodates single-port MAC.
- 6 m (20 ft) umbilical hose (longer available).
- Ballast for fresh or salt water (specify).
- Body dimensions – 178 mm diameter x 591 mm long (7 x 23.25 in.).
- Shipping weight – 20 kg (45 lb).

Model 7413, 7413F Detector Tripod



- For field measurements.
- Heavy duty construction.
- 7413 accommodates Slimline, MAC or Big MAC.
- 7413F accommodates Flanged, MAC or Big MAC.
- Detector element elevation (<30 cm to >120 cm).
- Built-in level.
- Shipping weight – 9 kg (20 lb).

